


Digital Removable Denture Prosthetics




01 >> Oral Scan

i-Vinci Intraoral Scanners



02 >> CAD Process

Choose Delux EPN tooth library in 3shape and exocad



03 >> Print Base

Denture Base Resin 3D Printer




05 << Great Result

Enjoy Perfect Smiles



04 << Teeth Bonding


Prefabricated Synthetic Polymer Teeth Digital Bond Kit



03 << Mill Base

Pink PMMA BLOCK Milling Machine

Implant-Supported Denture Prosthetics




01 >>

Take X-ray and decide the placement location for the dental implants



02 >>

Place implant inside the bone and add locator platform



03 >>

Take impressions and place supported denture onto model



05 <<

Final result



04 <<

Dentures with holes for implants



Functional Design

- > Universally adaptable and reliable for all concepts of occlusal patterns with multifunctional
- > Easy and intuitive with occlusion designed by the cogwheel principle
- > Efficient with precise and easy-to-find centric design
- > Semi-anatomical occlusal cusp design for 28° and 20°



Easy to Order

Web order: store.hugedentalusa.com

App order: TEETH ASAP

Contact sales: per business card info

Email and Phone:

info@hugedentalusa.com & 626-283-5808



DOWNLOAD
TEETH ASAP




Natural Aesthetics


- Multi-Layered perfection, designed for anterior excellence
- Harmonized layering of posterior teeth
- Smooth surface texture with aesthetics
- BL2, BL3, A00 and A0 available at request

Reliable Durability

- > Double cross-linked material, the Extended Polymer Network (EPN)
- > Higher molecular weight and wear-resistance
- > Lower water absorption and higher stain-resistance
- > More durable than conventional PMMA




PMMA Polymer



EPN Polymer

Performance of EPN Material^[1]



[1]: Test report from HUGE professional laboratory

Compatible with Digital Denture Workflow

- > Tooth library available in exocad & 3Shape
- > Enhanced workflow efficiency by 53%^[2]



[2]: Evaluation from HUGE's cooperating dental technicians

Mould Chart

Delux EPN

Let Your Smile Shine with NOBILDENT




UPPER ANTERIOR SQUARE




UPPER ANTERIOR OVOID




UPPER ANTERIOR TAPERING






UPPER ANTERIOR TAPERING





LOWER ANTERIORS




POSTERIORS 28°



POSTERIORS 20°



All measurements in milimeters



▲ SQUARE ■ TAPERING ● OVOID

Moulds


- > 11 upper anterior moulds
- > 7 lower anterior moulds
- > 12 posterior moulds

Shades

- > A1, A2, A3, B1, B2

Packing

- > Anterior: 6pcs/card
- > Posterior: 8pcs/card



Actual Size: All measurements in millimeters

SUGGESTED ARTICULATIONS

ANTERIORS		POSTERIORS	
UPPER ANTERIORS	LOWER ANTERIORS	POSTERIORS 28°	POSTERIORS 20°
S2	L2	30	30M
S3	L2	34	32M
S4	L3	34	32M
T2	L2	32	30M
T4	L4	32	30M
T6	L6	34	30M
T8	L8	34	32M
IO2	IL2/IL3	I33	32M
IS3	IL2/IL3	I33	32M
IT2	IL2/IL3	I33	32M
IT3	IL2/IL3	I33	32M